

Quality time for better care

"Quality time for better care" is the brand promise of Terumo Medical Care Solutions.



HOSPITAL CARE

AGANI™

Hypodermic Needle



Your needles of choice  
for a wide range of applications

# Smooth and easy puncture increase patient comfort

Make the administration of medication as painless as possible with the special geometry of the sharp lancet point and silicone-coated cannula that allows a smooth penetration and prevents any knocking effect<sup>1</sup>.



## Precisely shaped needle point with silicone coating.


Agani™ Hypodermic needles have a sharp lancet tip and an outside layer of silicone lubricant.

### The sharp point

- Facilitates the penetration of the skin
- Helps to minimize the trauma

### The silicone coating

- Provides good gliding properties
- Allows a smooth insertion



1 Data on file

# Life is easier with needles from Terumo

**Agani hypodermic needles provide you and your patients with the best possible levels of care and comfort.**

- Sharp lancet point – easy puncture
- Silicone coating – smooth needle movement
- Thin walled needles<sup>3</sup> – increased flow rate
- Short bevel options – reduced risk of vascular perforation
- Large portfolio – suitable for a wide range of applications



**Injection & Access**

**Contributing to Society  
through Healthcare**



**TERUMO EUROPE N.V.**

Terumo Corporation    Terumo Europe N.V.  
+81 3 6742 8500    +32 16 38 12 11

[www.terumo-europe.com](http://www.terumo-europe.com)

Your Terumo partner

All brand names are trademarks or registered trademarks owned by TERUMO CORPORATION, its affiliates, or unrelated third parties. Published by Terumo Europe N.V.

<sup>3</sup> Available on specific item codes

# Choose from a wide range from 18G up to 30G



## Agani™ Hypodermic Needles

The range includes a broad variety of needle lengths and diameters



**18G Cannula (ø 1.2 mm)**  
Length 38 mm  
Colour Pink



**19G Cannula (ø 1.1 mm)**  
Length 25, 38, 50 mm  
Colour Cream



**20G Cannula (ø 0.9 mm)**  
Length 25, 38 mm  
Colour Yellow



**21G Cannula (ø 0.8 mm)**  
Length 16, 25, 38, 50 mm  
Colour Green



**22G Cannula (ø 0.7 mm)**  
Length 32, 38 mm  
Colour Black



**23G Cannula (ø 0.6 mm)**  
Length 16, 25, 32 mm  
Colour Blue



**24G Cannula (ø 0.55 mm)**  
Length 25 mm  
Colour Purple



**25G Cannula (ø 0.5 mm)**  
Length 16, 25 mm  
Colour Orange



**26G Cannula (ø 0.45 mm)**  
Length 13, 23 mm  
Colour Brown



**27G Cannula (ø 0.4 mm)**  
Length 16, 19 mm  
Colour Grey



**30G Cannula (ø 0.3 mm)**  
Length 13 mm  
Colour Yellow



Does not contain  
components made of PVC  
or natural rubber latex.

For further ordering information, please refer to the Agani™ Safety Hypodermic needle product information card.

# Secure and easy handling increases safety and user comfort

The steeper point angle of the short bevel<sup>2</sup> shortens the length of the needle bevel, thereby reducing the risk of vascular perforation.

The outside diameter of the needle can be recognized by the ISO colour-coded hub.

- Easy identification of the required size
- Easy recognition of flashback through the semi-transparent hub

All needles come sterile, individually packed in colour-coded blister packaging

- Supplied in boxes of 100 needles
- Comprehensive information on the packaging




<sup>2</sup> Available on AN\*1838S1 and AN\*1938S1

The safety and  
comfort of your patients  
is in your hands



# Smooth and easy puncture increase patient comfort

Make the administration of medication as painless as possible with the special geometry of the sharp lancet point and silicone-coated cannula that allows a smooth penetration and prevents any knocking effect<sup>1</sup>.



## Precisely shaped needle point with silicone coating.

Agani™ Hypodermic needles have a sharp lancet tip and an outside layer of silicone lubricant.

- The sharp point**
  - Facilitates the penetration of the skin
  - Helps to minimize the trauma
- The silicone coating**
  - Provides good gliding properties
  - Allows a smooth insertion

<sup>1</sup> Data on file



Quality time for better care

<sup>2</sup> "Quality time for better care" is the brand promise of Terumo Medical Care Solutions.



# Life is easier with needles from Terumo

Agani hypodermic needles provide you and your patients with the best possible levels of care and comfort.

- Sharp lancet point – easy puncture
- Silicone coating – smooth needle movement
- Thin walled needles<sup>3</sup> – increased flow rate
- Short bevel options – reduced risk of vascular perforation
- Large portfolio – suitable for a wide range of applications



## Contributing to Society through Healthcare



TERUMO EUROPE N.V.

Terumo Corporation Terumo Europe N.V.  
+81 3 6742 8500 +32 16 38 12 11

www.terumo-europe.com

Your Terumo partner

All brand names are trademarks or registered trademarks owned by TERUMO CORPORATION, its affiliates, or unrelated third parties. Published by Terumo Europe N.V.

<sup>3</sup> Available on specific item codes

MP210GB-0919FK-1103.231E

## HOSPITAL CARE

# AGANI™

Hypodermic Needle



## Your needles of choice for a wide range of applications

The safety and comfort of your patients is in your hands



Secure and easy handling increases safety and user comfort

The steeper point angle of the short bevel<sup>2</sup> shortens the length of the needle bevel, thereby reducing the risk of vascular perforation.

The outside diameter of the needle can be recognized by the ISO colour-coded hub.

- Easy identification of the required size
- Easy recognition of flashback through the semi-transparent hub

All needles come sterile, individually packed in colour-coded blister packaging

- Supplied in boxes of 100 needles
- Comprehensive information on the packaging



<sup>2</sup> Available on AN\*1838S1 and AN\*1938S1

Choose from a wide range from 18G up to 30G



Agani™ Hypodermic Needles

The range includes a broad variety of needle lengths and diameters

|   |  |
|---|--|
| <br><b>18G Cannula (ø 1.2 mm)</b><br>Length 38 mm<br>Colour Pink         | <br><b>19G Cannula (ø 1.1 mm)</b><br>Length 25, 38, 50 mm<br>Colour Cream     |
| <br><b>20G Cannula (ø 0.9 mm)</b><br>Length 25, 38 mm<br>Colour Yellow   | <br><b>21G Cannula (ø 0.8 mm)</b><br>Length 16, 25, 38, 50 mm<br>Colour Green |
| <br><b>22G Cannula (ø 0.7 mm)</b><br>Length 32, 38 mm<br>Colour Black    | <br><b>23G Cannula (ø 0.6 mm)</b><br>Length 16, 25, 32 mm<br>Colour Blue      |
| <br><b>24G Cannula (ø 0.55 mm)</b><br>Length 25 mm<br>Colour Purple    | <br><b>25G Cannula (ø 0.5 mm)</b><br>Length 16, 25 mm<br>Colour Orange      |
| <br><b>26G Cannula (ø 0.45 mm)</b><br>Length 13, 23 mm<br>Colour Brown | <br><b>27G Cannula (ø 0.4 mm)</b><br>Length 16, 19 mm<br>Colour Grey        |
| <br><b>30G Cannula (ø 0.3 mm)</b><br>Length 13 mm<br>Colour Yellow     |  |

CE 0123 Does not contain components made of PVC or natural rubber latex.

For further ordering information, please refer to the Agani™ Safety Hypodermic needle product information card.