

Quality time for better care

"Quality time for better care" is the brand promise of Terumo Medical Care Solutions.



HOSPITAL CARE



Portfolio Guide

## Needles

### Agani™ needle


Sterile single use hypodermic needles with special geometry and sharp lancet point.

### Agani™ Safety hypodermic needle

Sterile single use safety engineered hypodermic needles to minimize the risk of an accidental needle stick.

### Dental needle

Sterile single use needle to be used with cartridge syringe suitable for dental anaesthesia and other purposes, available in a wide range of sizes.



## Syringes

### Terumo syringes

Sterile 3-part syringe for single use for the aspiration and/or injection of fluids for medical purposes.

### Syringes with needles

Sterile 3-part syringe for single use pre-connected with a sterile hypodermic needle.



## Miscellaneous

### Tourniquets

Tourniquets are an essential tool to prepare patient for venous access.




## Peripheral I.V. catheters

### Conventional I.V. Catheters

Short peripheral intravascular catheter with an ETFE<sup>1</sup>, FEP<sup>2</sup> or PUR<sup>3</sup> catheter tube and stainless steel introducer needle. Different models are available: straight, winged, winged and ported and with Surfflash™ system.

### I.V. catheters with integrated sharps protection

Short peripheral intravascular catheter with an introducer needle and integrated sharps protection technology.<sup>4</sup> Different models are available: straight, winged and winged and ported.




## Winged infusion sets

### Surflo™ winged infusion set

The Surflo winged infusion set (= Scalp Vein or SV-set) is used for intravenous injection.

### Surflo™ winged infusion set with protector

The Surflo winged infusion set with protector (= Scalp Vein or SV-set) is used for intravenous injection. The Surshield™ protector is intended to cover the exposed needle immediately after use to prevent needle stick injuries.



## Thermometer


Terumo professional thermometers guarantee two precise methods to measure body temperature and enable a long lifespan of approximately 10.000 measurements.



1 Surflo™ range  
2 Surflo™ S-Plus, -W and -WP range

3 Versatus™ and Surfflash™ range  
4 Surshield™ Versatus™ range

The safety and comfort  
of your patients is in your hands!



Peak performance  
even in little things - made easy!



**TERUMO EUROPE N.V.**

Terumo Corporation    Terumo Europe N.V.  
+81 3 6742 8500    +32 16 38 12 11

[www.terumo-europe.com](http://www.terumo-europe.com)

Your Terumo partner



All brand names are trademarks or registered trademarks owned by TERUMO CORPORATION, its affiliates, or unrelated third parties. Published by Terumo Europe N.V.

Please note that some products may not be available in all regions, check with your local Terumo contact.