

Terufusion[®]

Infusion System

Intensive Care Units



TERUMO

Terumo gives you a clear view on infusion

EASY HANDLING



All-in-one information at a glance

Healthcare professionals need to see at a glance which alarm or warning has been activated.

By colours and realistic illustrations, the Terfusion pump allows to differentiate instantly between high (red) and low (yellow) priority of alarms. Red alarms require an immediate intervention while yellow indications inform to prepare for an action in the (nearer) future.



Quick-Start of treatment

For healthcare professionals being able to start up or adapt the therapy quickly is not only an advantage but often a necessity for patient safety. The fast 'self check' of the Terufusion pump, combined with the simple interface and the unique rotary dial facilitates this.

If a patient condition deteriorates, it is necessary to adapt the flow rate quickly and simply. This is easily done by the turbo dial and by pressing the confirmation switch.



Risk reduction

On a daily basis healthcare professionals work with many different types of medication and dose regimes, which are quite complicated and might result in an error.

Using the Terufusion Drug Library reduces this risk. Drug types can be defined using the colour tag which will be shown on the screen.

Also the drug names can be programmed in Tall Man lettering to avoid making mistakes by choosing the correct drug.




Drug name with Tall Man letters	Can be confused with
DOBUtamine	DOPamine
buPROpion	busPIRone
chlorproMAZINE	chlorproPADIME
CEFOtaxim	CeftAZIdim
ATEnolol	AceBUTolol
cycloSPORINE	CYCLOserine
DimenhyDRINATE	DIPHENHYDRamine

Optimizing patient treatment

To optimize patient treatment a balanced and continuous infusion of multiple medications is imperative.

Internal communication between Terufusion SMART pumps stationed in a Terufusion communication rack system provides special functions like Piggy and Take Over Mode for a balanced continuous infusion.



Full focus infusion system

Terumo offers a comprehensive portfolio of solutions for the complete infusion line-up: user-friendly syringe and volumetric pumps and racks, as well as a wide range of highly accurate, easy-to-use disposables, such as an extended range of infusion sets and syringes.



Terumo Corporation
+81 3 3374 8111

Terumo Europe NV
+32 16 38 12 11

Your Terumo partner